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### Abstract

The new trends of contribution to various manufacturing technologies especially for the pharmaceutical market are introduced in terms of ‘Technology of Continuous Manufacturing with On-line Monitoring’ by continuous mixer based on FDA & ISPE. ‘Technology of dust explosion proof’ enables to prevent the disaster and risk for machine safety. ‘Technology of Isolator’ enables effective containment of hazardous and toxic materials and to avoid the scattering risk to the operators and the environment. ‘Technology of Freeze dryer’ and some compact machines would be useful for R&D and pharmaceutical examination.

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### Captions

- Fig. 1 Lift of major pharmaceutical company
- Fig. 2 Example of traditional tablet manufacturing process
- Fig. 3 Photo of modulomix
- Fig. 4 Example of continuous manufacturing with on-line monitoring
- Fig. 5 Example of two modulomix system
- Fig. 6 Modulomix (2units combined)
- Fig. 7 Overall efficiency
- Fig. 8 Isolator (Standard type)
- Fig. 9 Isolator (Compact type)
- Fig. 10 Isolator (Special type)
- Fig. 11 Example of isolator special-type
- Fig. 12 Hosokawa micron / bag in-bag out pulsjet collector BBP
- Fig. 13 Method of filter change
- Fig. 14 Flow of grinding equipment for cellulose derivative with explosion protection
- Fig. 15 Fine impact mill UPZ with pressure shock resistance (PSR 10bar)
- Fig. 16 Bag filter with PSR (Pressure shock resistance)
- Fig. 17 Active freeze dryer system
- Fig. 18 Evaluation of product by AFD
- Fig. 19 Image of AFD
- Fig. 20 Optisizer system
- Table 1 Forecast growth of pharmaceutical market